


# Iris Retractors / Hooks and Capsular Tension Rings



**For mechanical dilation of the pupil when pharmacological dilation is unsuccessful**

**Iris Retractor product include:**

- Stretches and retains the iris for maximum visibility during procedures such as phacoemulsification, retina and vitreous surgery.
- Easy to use, safe, fast alternative to reduce the surgical time.
- Suitable for aphakic, pseudophakic and phakic eyes.
- Ideal temporary implant.

**Product code:** IH100 (box of 5)

**Enhance safety and effectiveness during phacoemulsification and posterior chamber intraocular lens (IOL) implantation**

Stabilise weakened, broken, or missing zonules during cataract extraction using phacoemulsification and continuous curvilinear capsulorhexis techniques in the human eye.



**CTR product benefits:**

- SCircular expansion and stabilisation of the periphery of the capsular bag.
- Prevents vitreous herniation from weak zonular area.
- Safe IOL centration in eyes with zonular dehiscence.
- Reduced risk of capsular fibrosis and shrinking.

Model	CTR 11	CTR 12	CTR 13
Overall Diameter	11.0 mm	12.0 mm	13.0 mm
Insertion Diameter	9.7 mm	10.4 mm	11.4 mm
Product code	CTR11 W CTR11 B	CTR12 W CTR12 B	CTR13 W CTR13 B

CTR Range available in Clear (W)  
or Blue (B) PMMA CQ-UV

Available in the UK from



Tel +44 (0)161 776 7627 or email [sales@surgitrac.co.uk](mailto:sales@surgitrac.co.uk)  
Visit us online at [www.surgitrac.com](http://www.surgitrac.com)

10 Wharfside Business Park, Irlam Wharf Road, Irlam,  
Manchester, M44 5PN UK

