

Bio Blue™ 90



Acid Blue 90 Ophthalmic Solution

Bio Blue™ is a sterile solution of Trypan Blue, a biocompatible intraocular solution used in the human eye during ophthalmic surgery.

Bio Blue™ 90 is a sterile, non-pyrogenic, isotonic solution of Acid Blue 90 (a non-fluorescent Triphenylmethane) stain for use in the human eye.

- Ready to use.
- No mixing or filtration is required.
- Easy to use and remove.



Indications include:

- An adjunct to vitreoretinal surgery
- Enhances visualisation of the Epiretinal Membrane during Macular Hole and Pucker surgery.
- Enhances visualisation and facilitates ILM peeling.
- Capsular stain for Continuous Curvilinear Capsulorhexis.

Specifications:

- Bio-Blue 90 Acid Blue, buffered aqueous base.
- 0.5 ml BD luer-lock glass syringe.
- 27 G CE marked cannula.
- Ambient storage – no refrigeration required.
- Non pyrogenic and latex free.

Supplied in 0.5 ml PFS.

Product code: Bioblue90 (single unit)

Available in the UK from

Surgitrac® Instruments UK Limited

10 Wharfside Business Park, Irlam Wharf Road, Irlam, Manchester, M44 5PN UK

Tel: +44 (0)161 776 7626 Fax: +44 (0)161 776 7616

sales@surgitrac.com www.surgitrac.com

