

Trypan Blue Ophthalmic Solution 0.06 % W/V

Bio Blue™ is a sterile solution of Trypan Blue, a biocompatible intraocular solution used in the human eye during ophthalmic surgery.

Selectively stains connective tissue structures such as, the anterior lens capsule of the crystalline lens.

Trypan blue is intended to be applied directly on the anterior lens capsule - staining any portion of the capsule which comes in contact with the dye. The dye does not penetrate the capsule, permitting visualisation of the anterior capsule in contrast to the non-stained lens cortex and inner lens material.



**For
Cataract
surgery**



Indications include:

- Phacoemulsification.
- Small Incision Cataract surgery.
- Extra capsular Cataract surgery.
- Anterior / Posterior Capsulotomy.
- DMEK and DALK.

Specifications:

- Bio Blue™ Trypan Blue Ophthalmic Solution 0.06% W/V.
- 1.0 ml sterile pre-filled BD luer-lock glass syringe.
- Sterile 27G CE marked cannula.
- Ambient storage – no refrigeration required.
- Non pyrogenic and latex free.

Product code: Bioblue (single unit)

Available in the UK from

Surgitrac® Instruments UK Limited

10 Wharfside Business Park, Irlam Wharf Road, Irlam, Manchester, M44 5PN UK

Tel: +44 (0)161 776 7626 Fax: +44 (0)161 776 7616

sales@surgitrac.com www.surgitrac.com

