

Ophthalmic Cannulas



- Single-use
- Available sterile or non-sterile
- Can be included in procedure packs



Hub Colour Guide

Colour	Gauge	
Natural* (Clear)	30G	
Grey	27G	
Brown	26G	
Orange	25G	
Purple	24G	
Blue	23G	
Black	22G	
Green	21G	
Yellow	20G	
Cream	19G	

Anaesthesia Cannulas

Retro bulbar Atkinson

Administer topical anaesthetics into the muscle cone of the eye.

Retro bulbar Atkinson 23G, 25G, 27G

Retro bulbar non -Atkinson

for depositing anaesthesia outside the muscle cone around the globe.

Retro bulbar non -Atkinson 25G, 26G, 27G

Peribulbars Atkinson 23G, 25G, 27G

Peribulbars non -Atkinson 23G, 25G, 27G, 30G

Rous sub-tenon / Para bulbar needle -

For administering anaesthetic agent in the posterior sub-Tenon's space.

Rous sub-Tenon / Para bulbar needle 19G, 21G



Anterior chamber cannulas

Regular wall (angled, curved, straight) 19G, 20G, 21G, 23G, 25G, 27G, 30G


Thin wall (angled, curved, straight) 21G, 25G



Capsule polisher

Polishing and cleaning residual cortex in the capsular bag before IOL insert.

Sandblasted 21G, 23G
Silicone tip 23G, 27G



Capsulorhexis Cystomes

For Capsulorhexis sharp tipped Cystomes.

90 degree Tip is angled at 90 degree then rotated 90

For Capsulorhexis sharp tipped Cystomes



Hydro dissection and Hydrodelineation

Designed to deliver fluids to facilitate separation of the cortex from the nucleus.

Angled flat bend 19G, 23G, 25G, 27G, 30G

Cortical cleaving 27G


Curved 27G

Rotator 27G

J-Shaped 25G, 27G, 30G

Kellen-Curved 25G, 27G

Chang 27G



Irrigating Cystomes

For performing Capsulorhexis.

Angled 25G, 27G

Curved 25G, 27G

Formed 23G, 25G, 27G, 30G

Straight 30G



Contact your local distributor for more information, or a quote.

Surgitrac® Instruments

Tel: +44 (0)161 776 7626

For UK enquiries contact info@surgitrac.co.uk Non UK customers contact info@surgitrac.com

www.surgitrac.com



0126

