










Instructions for use (IFU)

Surgitrac Instruments sterile single-use instruments

Device name	Single-use ophthalmic surgical instruments
Device codes	<p>SC01,SC05,SC07,SC08,SC09,SC10,SC11,SC11P, SC12,SC13,SC14,SC15,SC16,SC18,SC19,SC20,SC21 SC22,SC23,SC24,SC25,SC29,SC30,SC31,SC32,SC33,SC34, SC35,SC36,SC37,SC38,SC39,SC42, SC43,SC44,SC45,SC45 SC46,SC47,SC48,SC49,SC50,SC51,SC52,SC53,SC54,SC55 SC56,SC58,SC59,SC60, SC61, SC62, SC63, SC64, SC65, SC71 SC72, SC73, SC80, SC81, SC85, SC90, SC91, SC92, SC93 SC140, SC141.</p> <p>SR10,SR11,SR21,SR31,SR32,SR34,SR40,SR41,SR43,SR51,SR60 SR61,SR71,SR90,SR91.</p> <p>SG01,SG02,SG03,SG03A,SG04,SG04P,SG04S,SG05,SG06,SG09 SG10,SG11,SG12,SG13,SG14,SG15,SG17,SG19,SG20,SG22,SG23 SG24,SG92,SG93</p> <p>CR01,CR02,CR03,CR05 CS01,CS02,CS05,CS05, SA11,SA12,SA13</p> <p>PD01,PD02,PD03 CA03,CA04</p> <p>RMPL, RMPL2,RMPS,RMPS2</p>
Item specifics	For use in ophthalmic surgery
Intended use	For use in the human eye.
For use by	Only for use by suitably trained medical professionals.
Warning	Used surgical instruments are contaminated. Disease may occur from injury. Dispose of in a puncture resistant container. Refer to local regulatory guidelines
Risks, contradictions and side effects	Incorrect use may cause blindness.
	<p>Manufacturer Sihong Surgitrac Co. Ltd No. 2, Xihu Road, Sihong Industrial Park, Jiangsu, 223900 www.surgitrac.com info@surgitrac.com +44 (0) 161 775 7626</p>
	<p>Single-use only. Do not re-use.</p>
	Sterilized using radiation.

	Keep dry.
	Keep away from Sunlight.
  eIFU	Consult instructions for use Electronic version available at www.surgitrac.com
	Caution see instructions for use.
	Do not use if package is damaged.